



Product Code	27708
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184716
CND	Y122106
EAN/UPC	8023279277081
UMDNS	16949
PARAF	987471493

Description

PAEDIATRIC WHEELCHAIR - 35 cm seat - black tissue

- single cross bar
- padded flip-up armrest
- fixed footrest, flip-up foot plate with leg strap
- 15 cm (6") solid front castors
- 56 cm (22") pneumatic rear wheels
- locking brakes on rear wheels
- folding back
- large pocket on the back

Removable padding with velcro fixing

Technical Specifications

- Total width: 52 cm
- Folded dimensions: 32 cm
- Seat width: 35 cm
- Rear wheel diameter: 56 cm
- Front wheel diameter: 15 cm
- Seat height: 48 cm
- Total wheelchair height: 90 cm
- Total wheelchair length: 97 cm
- Seat depth: 36 cm
- Backrest height: 40 cm
- Max load: 75 kg
- Kit size: 85 x 33 x 73 cm
- Net Weight: 13 Kg
- Gross Weight: 15 Kg