

ADJUSTABLE SWIVEL STOOL



Product Code	27733
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279277333
PARAF	977812078

Description

ADJUSTABLE SWIVEL STOOL WITH HANDLE

- swivel seat with 90 degree position lock at every right angle location
- non-slip central seat pad and soft under-covers
- load 135 kg (300 lbs)

Adjustable seat 54.4 - 71.6 cm with 8 height adjustments, by 2.5 cm (1") increment