



Product Code	27738
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

CND	Y091215
EAN/UPC	3701127751116
GMDN	91218

Description

CLIPUP® TOILET RAISED - autoclavable

ClipUp® can be installed safely on any type of toilet bowl.
Innovating fixation system with ratchet.
Design with front cutout provides optimal comfort and facilitates personal hygiene.
Max Load 150 kg.

Delivered with 6 antislip SBS pads.
Patented fastening system with ratchet.

Made in France.

Technical Specifications

Size:

- width 34 to 37 cm
- depth 38.5 cm
- height 14 cm

Packing: 45.5x16.5xh 39.5 cm
Seat opening (depth x width): 26x16 to 19 cm
Toilet elevation: 10 cm
Weight: 0.872 kg
Material: 100% polypropylene