

RAISED TOILET SEAT with reclining armrest - height 12,5 cm



Product Code	27751
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2566853
CND	Y091215
EAN/UPC	8023279277517
PARAF	976329728
GMDN	31059

Description

RAISED TOILET SEAT height 12.5 cm

With fixing system (back and front screw) and anodized aluminium reclining armrest.

Max load 130 kg (286.6 lbs)

Packing: 410 x 165 x 468 mm

Instructions: GB, IT, FR, ES, DE, GR

Made in Italy.