

# RAISED TOILET SEAT with reclining armrest - height 12,5 cm



Product Code	27751
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2566853
CND	Y091215
EAN/UPC	8023279277517
PARAF	976329728
GTIN	08023279277517
GMDN	31059

## Description

RAISED TOILET SEAT height 12.5 cm

With fixing system (back and front screw) and anodized aluminium reclining armrest.

Max load 130 kg (286.6 lbs)

Packing: 410 x 165 x 468 mm

Instructions: GB, IT, FR, ES, DE, GR

Made in Italy.