

CLIPPER RAISED TOILET SEAT - height 11 cm



Product Code	27756
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2044555
CND	Y091203
EAN/UPC	3401077204822
PARAF	976329716

Description

CLIPPER RAISED TOILET SEAT - 11 cm

Easily adapt to any toilet bowl.

Totally made of injected plastic.
4 fixation clips to secure strongly the seat on the toilet bowl.
Max load 185 kg (407 lb).

Made in France.