



Product Code	27800
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831159
CND	V0899
EAN/UPC	8023279278002
PARAF	903113456
GMDN	35892

## Description

### PATIENT TROLLEY

Made of chrome plated steel tube. Comes complete with mattress made of washable self extinguishing plastic material.  
 Mounted on 4 castors Ø 100 mm two of which with brakes.  
 Size: 180x54xh 80 cm.  
 Load: 110 kg.

Adjustable backrest up to a 60°

Multilingual instructions: GB, FR, IT, ES

Made in Italy.