

# PATIENT TROLLEY - REMOVABLE TOP



|               |                |
|---------------|----------------|
| Product Code  | 27801          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | No  |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2831159        |
| CND        | V0899          |
| EAN/UPC    | 8023279278019  |
| PARAF      | 903113468      |
| GTIN       | 08023279278019 |
| GMDN       | 35892          |

## Description

PATIENT TROLLEY - removable top

Made of chrome plated steel tube. Comes complete with mattress made of washable self extinguishing plastic material.  
Mounted on 4 castors Ø 100 mm two of which with brakes.  
Size: 180x54xh 80 cm.  
Load: 110 kg.

Adjustable backrest up to a 60°

Multilingual instructions: GB, FR, IT, ES

Made in Italy.