

PROFESSIONAL PATIENT TROLLEY with side rails and oxygen cylinder holder



Product Code	27804
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2263972
CND	V0899
EAN/UPC	8055136920250

Description

PROFESSIONAL PATIENT TROLLEY with adjustable head section, side rails and oxygen cylinder holder

Innovative, stylish comfortable trolley easy to use.
Trolley made of sturdy tubular frame for optimal stability.
Supplied with 4 wheels Ø 200 mm (2 with brake), bumpers, push handle.
Detachable frame with adjustable head section.
Size: 196x76xh 75 cm.
Net size: 193x68 cm
Maximum load: 200 kg.
Weight: 70 kg.

Made in Italy.