

PROFESSIONAL RX PATIENT TROLLEY with side rails and oxygen cylinder holder



Product Code	27805
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2264002
CND	V0899
EAN/UPC	8055136920243

Description

PROFESSIONAL PATIENT TROLLEY - RADIOTRSPARENT FRAME with adjustable head section, side rails and oxygen cylinder holder

RADIOTRSPARENT FRAME

Position of vertical columns and radiotransparent frame allow the structure to have a wide radiotransparent surface.

Innovative, stylish comfortable trolley easy to use.

Trolley made of sturdy tubular frame for optimal stability.

Supplied with 4 wheels Ø 200 mm (2 with brake), bumpers, push handle.

Fixed frame, 12 mm thickness, with adjustable head section.

Size: 196x76xh 75 cm.

Net size: 193x68 cm

Maximum load: 200 kg.

Weight: 70 kg.

Made in Italy.