



Product Code	27826
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	341574
CND	V0899
EAN/UPC	8023279278262

## Description

Infusion stand

Made of chrome plated steel.  
Four chrome plated steel hooks.  
Five spoke stainless steel base Ø 60 cm with castors Ø 50 mm, two with brake.  
In conformity with MDR.

Made in Italy.

## Technical Specifications

Adjustable height (cm): 160-250  
Load (kg): 30 (5 kg each hook, 10 kg for infusion pump frame)  
Weight (kg): 10  
Infusion pump frame (cm): 39 x 75  
Packaging size (cm): 185 x 18 x 43