

SHOWER MATTRESS



Product Code	27836
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279278361

Description

SHOWER MATTRESS

It allows bed-ridden, elderly or disabled patients to be washed directly on their bed, without any need for transfer/lift. The patient can be showered with water coming from the tap via an extension hose or using the special pressurized water tank provided with flexible hose and shower head. Once finished, the water is drained into a suitable container. Complete with plug and flexible hose 6.4 m for quick water drain.

Made in Italy.