

I.V.STAND ON 5 WHEELS TROLLEY - stainless steel - 4 s/s hooks



Product Code	27839
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2753826
CND	V9099
EAN/UPC	8023279278392
PARAF	977812142
GMDN	36069

Description

Stainless steel AISI 304 IV stand

Adjustable height: 160-210 cm

Weight: 3.3 kg

Base material: Stainless steel

Hook (number): 4

Hook material: Stainless steel

Load kg Per hook: 1.2 kg

Total load: 4.8 kg

Castors with brake: 1

Castors Antistatic: -