

I.V.STAND ON 5 WHEELS TROLLEY - stainless steel - 4 plastic hooks



Product Code	27844
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2823532
CND	Z129099
EAN/UPC	8023279278446
PARAF	977812179
GTIN	08023279278446
GMDN	36069

Description

Stainless steel AISI 304 IV stand

Adjustable height: 135-215 cm

Weight: 2.9 kg

Base material: Stainless steel

Hook (number): 4

Hook material: Plastic

Load kg Per hook: 1.2 kg

Total load: 4.8 kg

Castors with brake: 1

Castors Antistatic: 5