



Product Code	27857
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2831153
CND	V0899
GMDN	36069

## Description

### HYPODERMOCLYSIS STAND

Stand with hydraulic height adjustment for hypodermoclysis, for washing and surgical irrigation. Hydraulic pump with pedal for height adjustment and gripping handle for easy transport.

LOAD 28 KG

4 stainless steel AISI 304 hooks, with antirelease safety device

Stainless steel adjustable gripping handle for easy transport

Stainless steel AISI 304 stand Ø 60 mm

Pedal driven shock absorbing pump.

Height adjustable: 147-247 cm

Weight: 22 kg

Stainless steel base with 5 pivoting castors Ø 75 mm (2 antistatic with brake, 3 standard) and grey vipla bumpers.

Made in Italy.