

# SUPER EASY CART with shelf 32x26 cm



Product Code	27870
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	173908
CND	V0880
EAN/UPC	8023279278705
GMDN	173908

## Description

SUPER EASY CARTS - with handle

- Cart + shelf + base
- Shelf size (cm): 32x26

A new stylish and quality cart you can build according to your needs. Carts with base and 1 or 2 sturdy shelves, aluminium columns and 4 twinned castors Ø 65 mm (2 with brake). Illustrated assembly instructions.

BUILD UP YOUR CART, YOU CAN ADD: EXTRA SHELF, BASKET, INFUSION STAND, CARDIO ARM.

Made in Italy.

## Technical Specifications

Dimensions: 47x42x79 cm  
Load capacity: total 120 kg - shelf 10 kg  
4 twinned castors Ø 65 mm (2 with brake)  
Material:  
- extrusions in aluminium  
- plastic parts in ABS  
Shelf: steel sheet epoxy powder coated