

EASY CART - 2 shelves 40x36 cm + base



Product Code	27880
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	173918
CND	V0880
EAN/UPC	8023279278804
GMDN	33398

Description

EASY CARTS

- 2 shelves 40x36 cm + base
- Height: 75 cm

A complete line of stylish and sturdy carts with 2, 3 or 4 shelves. Adjustable height of sturdy shelves, aluminium columns and 4 twinned castors Ø 65 mm (2 with brake). All shelves with anti-fall edges. Designed for medical devices. Illustrated assembly instructions.

Made in Italy.

Technical Specifications

Dimensions (WxDxH): 50x45x75 cm

Load capacity:
- total load 140 kg
- shelf 25 kg

Shelves: 40x36 cm

Castors: 4 twinned castors Ø 65 mm (2 with brake)

Material:
- vertical aluminium extrusions
- plastic parts in ABS

Shelf: steel sheet epoxy powder coated