

# PRO CART - 2 shelves, keyboard, monitor support



Product Code	27887
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	109277
CND	V0880
EAN/UPC	8023279278873

## Description

PRO CART 2 shelves + keyboard and monitor support

All professional carts are supplied with cable pipe inside lateral columns.

Made in Italy.

## Technical Specifications

Load capacity: - total load 240kg - shelf 35 kg

Base: 60.8x50 cm

Shelves: 2 shelves 50x50 cm + base

Material:

- vertical extrusions in aluminium

- plastic parts in ABS