

# SMART CART - 2 shelves 50x42 cm + base



Product Code	27889
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2126913
EAN/UPC	8053736710585

## Description

SMART CART - 2 shelves 50x42 cm + base

Sturdy carts with painted powder steel base and aluminium frame. 2 shelves with protection edge on back, adjustable at different height.

Suitable for all Gima electrosurgical units.

Illustrated assembly instructions.

First shelf can be mounted flat or tilted  $\pm 10^\circ$

4 twinned pivoting castors (2 with brake)  $\varnothing$  65 mm made of thermoplastic rubber and ABS.

Made in Italy.

## Technical Specifications

Total load: 80 kg  
Load x shelf: 15 kg  
Size mm: 582x500xh 780  
Weight: 15 kg