



Product Code	27896
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	173938
CND	V0880
EAN/UPC	8023279278965

## Description

SMART CARTS - grey/white RAL 7035

- 3 shelves 40x36 cm + base
- Height: 88 cm

Carts with base and 2, 3, or 4 wide, sturdy shelves in 3 different widths 40, 50 or 60 cm.  
Aluminium columns and 4 twinned castors Ø 65 mm (2 with brake).  
First shelf can be mounted flat or tilted  $\pm 10^\circ$ .  
Designed for medical devices and versatile for other uses.  
Illustrated assembly instructions.

All shelves have protection edge on back.  
Shelves adjustable at different height.  
Shelves made with coat of paint sheets.

Made in Italy.

## Technical Specifications

Dimensions (cm): 48,2x50x88

Load capacity: total load 80 kg, shelf 15 kg

Material:  
- vertical aluminium extrusions  
- plastic parts in ABS

Shelf: steel sheet epoxy powder coated

Total weight: 15 kg