

VALUE CABINET - 1 door - tempered glass



Product Code	27898
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279278989
---------	---------------

Description

VALUE CABINET - 1 door

Cabinet made of enamelled steel sheet with tempered glass door.

Hinged glass door with lock and key.

Three stainless steel shelves adjustable in height.

Delivered in kit form.

Size: 52x35x138 cm

Pack size: 65x26x134 cm

Weight: 31 kg

Made in Italy.