

# VALUE CABINET - 2 doors - tempered glass



Product Code 27899

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279278996

## Description

### VALUE CABINET - 2 doors

Cabinet made of enameled steel sheet with tempered glass doors.

Hinged glass door with lock and key.

Three stainless steel shelves adjustable in height.

Delivered in kit form.

Size: 75x37x164 cm

Pack size: 72x24x156 cm

Weight: 45 kg

Made in Italy.