

# VALUE CABINET - 4 hinged doors - tempered glass



Product Code 27900

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279279009

## Description

### VALUE CABINET - 4 doors (hinged)

Cabinets made of enameled steel sheet with 2 overlapped compartments.

Upper compartment with 2 tempered glass hinged doors, lock with key and 2 stainless steel, height adjustable shelves. Inferior compartment with 2 hinged doors made of enameled steel plate and 1 enameled steel sheet shelf.  
Delivered in kit form.

Size: 75x37x164 cm

Pack size: 72x24x156 cm

Weight: 50 kg

Made in Italy.