

VALUE CABINET - 4 doors - tempered glass



Product Code 27902

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279279023

Description

VALUE CABINET - 4 doors (2 hinged, 2 sliding)

Cabinets made of enameled steel sheet with 2 overlapped compartments.

Upper compartment with 2 tempered glass sliding doors, lock with key and 2 stainless steel, height adjustable shelves. Inferior compartment with 2 hinged doors made of enameled steel plate and 1 enameled steel sheet shelf.

Delivered in kit form.

Size: 79 x 36 x 186 cm

Pack size: 88 x 36 x 96 cm

Weight: 56 kg

Made in Italy.