



Product Code	27911
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	Not required
--------	--------------

EAN/UPC	8023279279115
---------	---------------

PARAF	903608305
-------	-----------

## Description

Hanging cabinet to keep medicine with 2 glass tempered doors and 2 tempered glass shelves adjustable on master rack to be hanged by two holes on the back. Supplied with lock.

Made in Italy.

## Technical Specifications

- Doors: 2
- Adjustable shelves: 2
- Weight: 16 kg
- Size: 60 x 30 x h 60 cm