



Product Code	27912
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	Not required
--------	--------------

EAN/UPC	8023279279122
---------	---------------

GTIN	08023279279122
------	----------------

Description

Hanging cabinet with 2 lockable doors.
Inside division in 2 compartments left side with 2 removable shelves, right side with 1 removable shelf.

Made in Italy.

Technical Specifications

- Doors: 2
- Adjustable shelves: 3
- Weight: 18 kg
- Size: 80 x 33 x h 45 cm