

TRUSTEEL NON ABSORB. SUTURE gauge 4 circle 1/2
needle 48 mm - 45 cm - natural



Product Code	22792
Unit of sale	box of 48 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	2381046
CND	H0102010199
EAN/UPC	8023279227925

Description

TRUSTEEL - non absorbable monofilament steel surgical suture

Equivalent Ethicon brand: ETHISTEEL

- Gauge: 4
- Metric: 6
- Colour: Natural
- Circle: 1/2
- Needle profile: Conventional cutting
- Needle length: 48 mm
- Suture length: 45 cm

Technical details:

- monofilament stainless steel LVM 316 grade wire
- optimal compatibility with body tissue
- low or no tissue reactivity
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- orthopedic surgeries
- sternum closure

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.