



| | |
|---------------|----------------|
| Product Code | 27921 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2255198 |
| CND | V0899 |
| EAN/UPC | 8055136920441 |

Description

Structure in oval tube, painted with epoxy powder metal grey. White shelves with PVC edge, guard-rail on 3 sides.

Size 50x60xh80 cm.
Max load: 30 kg.