



Product Code	27935
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279279351
---------	---------------

Description

ISO SERVICE TROLLEY - empty

Service trolley for "ISO" trays or baskets (not included) made of epoxy powder painted steel. Supplied with 5 height adjustable runners for ISO trays/baskets.

Size: 65x47xh 100 cm

Weight: 12.5 kg

Pivoting castors Ø 100 mm (2 with brake)

Made in Italy.