



Product Code	27941
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279279412

## Description

COMBO MEDICINE CABINETS 96x66xh 195 cm

With 30 medicine compartments one side and 6 ISO trays 1 shelf and 3 drawers other side

- made of bilaminated board 20 mm thick
- epoxy powder painted steel base frame with 4 adjustable feet Ø 50 mm

### Right side

- 1 external lockable door turning through 180° with 10 medicine compartments each
- perimetrical anti-dust gaskets
- inside 1 container with 10 medicine compartments each side, turning through 90°, provided with lock

### Left side

- glass door, internal side, with runner for ISO trays
- 3 drawers on the bottom

Supplied with 6 ISO trays different heights (2xh 5 cm, 3xh 10 cm and 1xh 20 cm) and 1 shelf.

Made in Italy.