

# MULTIFUNCTION STORAGE CABINET - light grey



Product Code	27959
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC 8023279279597

## Description

### MULTIFUNCTION STORAGE CABINET - light grey

Multifunction storage cabinet has a structure in powdercoated aluminium and profiles with resistant and durable melamine-faced wood particle panels.

- 2 doors with powder-coated aluminium edges with 5 mm tempered glass that rotate on 180° hinges with lock
- 4 adjustable shelves made of tempered glass

At the base 4 powder coated steel feet with adjustable plastic insoles, 200 mm height to facilitate cleaning.

Size: 750x400xh 1700 mm