



Product Code	27961
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279279610

## Description

### SHELL PLUS MEDICINE CABINET - white

Structure made of powder-coated aluminium with melamine panels.

- 2 doors made of 20 mm thick wood particle panels edged in white ABS, revolving on 270°, hinges complete with glass-filled nylon lock and handles
- 20 unbreakable and recyclable transparent compartment trays (10 on each door)
- 5 adjustable shelves with a narcotic cabinet with lock
- 2 revolving panels with 40 unbreakable and recyclable transparent trays (20 on each panel).

4 powder coated steel feet Ø 60 mm with adjustable plastic insoles.  
Size: 930x600xh 1980 mm