

TRANSPARENT PLASTIC ISO DRAWER 600x400x100 mm - open



| | |
|--------------|-------|
| Product Code | 27965 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

| | |
|------|-------------------|
| Type | No medical device |
|------|-------------------|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|--------|--------------|
| UK-REP | Not required |
|--------|--------------|

| | |
|--------|--------------|
| CH-REP | Not required |
|--------|--------------|

| | |
|---------|---------------|
| EAN/UPC | 8023279279658 |
|---------|---------------|

Description

Plastic tray for ISO trolleys (codes 27934-35-36).

Intelligent solutions for storing all kind of drugs, bandages, accessories and small components in a space saving and clear manner.

- Type: open
- Dimensions: 600 x 400 x 100 mm
- Material: Transparent PC