



Product Code	27990
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279279900
PARAF	903608560

Description

<div>STUDIO LIGHT - 100 W</div><div>
 </div><div>Halogen lamp with aluminium tubular frame and flexible chromed steel arm. Plastic base Ø 60 cm with plastic wheels Ø 50mm (1 with brake). 220-240 V, 50-60 Hz.</div><div>
 </div><div>Bulb (connection E27) not included.</div><div>
 </div><div>Flexible arm 38 cm</div><div>Aluminium dome Ø 15 cm</div><div>Net weight: 3 kg</div><div>Packing size: 115x17x57 cm</div><div>Fixed height 100 cm</div><div>
 </div><div>Made in Italy.</div><div>
 </div><div>
 </div>