



Product Code	28008
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279280081
---------	---------------

Description

AMIRA chair - 2 motors - white

Electrically controlled chair to adjust separately height and backrest.
Manual adjustment of angle of the seat inclination and footrest.
Supplied with cushion.

Face hole behind cushion

Size: 185 x 69 x h 64-86 cm

Net weight: 62 kg

Max load: 175 kg

Operating power: 100-240 V - 50/60 Hz.

Angle adjustment of the seat inclination 20° (manual)

Footrest angle adjustment 0° to 85° (manual)

Footrest manual extension

Height adjustment 64-86 cm