

AMIRA CHAIR - electric 2 engines - blue



Product Code	28009
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280098

Description

AMIRA chair - 2 motors - blue

Electrically controlled chair to adjust separately height and backrest.
Manual adjustment of angle of the seat inclination and footrest.
Supplied with cushion.

Face hole behind cushion
Size: 185 x 69 x h 64-86 cm
Net weight: 62 kg
Max load: 175 kg
Operating power: 100-240 V - 50/60 Hz.
Angle adjustment of the seat inclination 20° (manual)
Footrest angle adjustment 0° to 85° (manual)
Footrest manual extension
Height adjustment 64-86 cm