

NEW PODOLOGY MECHANICAL CHAIR - white



Product Code	28021
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279280210
---------	---------------

Description

This chair can move up and down and rotate.
Backrest and footrests are adjustable.
Heavy round base provides maximum stability.

Size: 160 x 73.5 x h 51-65 cm

Net weight: 41 kg

Max load: 200 kg