

# PODOLOGY ELECTRIC CHAIR 3 motors - blue



Product Code	28025
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280258

## Description

### ELECTRIC PODOLOGY CHAIR - 3 motors - blue

Electrically controlled chair to independently adjust height, backrest and seat inclination, by foot switch or push button panel. Head rest can be inclined, arm rests are movable.

Size: 186 x 86 x h 58-89 cm

Net/Gross weight: 72/92 kg

Max load: 175 kg

Operating power: 220-240 V - 50/60 Hz.