

TRULON NON ABSORB. SUTURE gauge 3/0 circle 3/8
needle 24 mm - 45 cm - blue



Product Code	22804
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2920515
CND	H010201010102
EAN/UPC	8023279228045
PARAF	952040160
GMDN	16584

Description

TRULON - polyamide nylon non absorbable monofilament suture
Equivalent Ethicon brand: ETHILON / NUROLON

- Gauge: 3/0
- Metric: 2.0
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 24 mm
- Suture length: 45 cm

Technical details:

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.