

NEFERTITI ELECTRIC CHAIR 3 motors - white



Product Code	28046
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280463

Description

EFERTITI ELECTRIC CHAIR - 3 motors - white

Electrical chair with very stable structure and high comfort.
Specially designed to optimize the workspace due to 240° rotation.
3 motors control height, backrest and seat inclination by foot switch or push button panel.
Footrest adjustable by gas lift pump, reversible arms and adjustable headrest.
Operating power: 220-240 V - 50/60 Hz.
Manual: GB, FR, IT, ES.

Size: 186x86xh 58-88 cm
Net weight: 79 kg
Max load: 175 kg
58-89 cm Height adjustment