



Product Code	28047
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280470

Description

SABA ELECTRIC chair - 4 motors - white

Electrically controlled 4 motors chair to adjust separately angle of the backrest, seat inclination, footrest and lifting of the whole chair.
 Hand controller and built-in foot controller which allow operator to adjust chair.
 High-grade upholstery leather easy to clean and fresh feeling.
 Headrest with face hole.
 Supplied with removable cushion and remote controller.
 Metal base for high stability.
 Operating power: 100-240V - 50/60Hz
 Manual: GB, FR, IT, ES.

Size: 185 x 88 x h 63-90 cm
 Net weight: 98 kg
 Max load: 175 kg
 Height adjustment of the whole chair 63-90 cm
 Angle adjustment of the seat inclination 0° to 20°
 Angle adjustment of the backrest 0° to 80°
 Footrest angle adjustment 0° to 65°