



Product Code	28051
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2741554
CND	Z120105
EAN/UPC	6933257923423

Description

AROS BATTERY SUCTION UNIT - 1 I

This lightweight and easy to carry suction unit is ideal for home and emergency use. Modern and portable design, it's equipped with a rechargeable lithium battery that delivers up to 60 minutes of use on a full charge. The suction unit features a dual jar hook system, allowing quick changeovers without interrupting the procedure. It also supports disposable liner bags (optional) to simplify cleanup and minimize the risk of cross-contamination.

- upgraded vacuum gauge and adjustment knob for accurate suction control

- built-in overflow protection system safeguards the internal pump

- real time battery indicator to monitor power levels

Technical Specifications

Air flow: 27 ±4 l/min

Operating cycle: continuous use with AC or 60 min with fully charged battery

Vacuum regulating range: 20kPa to maximum vacuum

Max vacuum: 85 ±5 kPa

Noise level: ≤ 65 dB(A)

Power supply: rechargeable lithium battery, 2600 mAh

Charging time: 4 hours

Operating voltage: 100-240 V, 50/60 Hz

Input power: 150 VA

Size: 480x189x285 mm

Weight: 4 kg

Protection: IP21, Class II

Standard accessories

Autoclavable jar 1,000 ml

Hose

Filters

Suction catheter

Fuse

Power cord

Suction tube 2m

User manual GB, FR, IT, ES, PT, DE NL, PL, SE, HU, RO, GR