



| | |
|---------------|----------------|
| Product Code | 28060 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|------------|---------------|
| RDM (NSIS) | 1734414 |
| CND | Z12159002 |
| EAN/UPC | 8054619104644 |
| PARAF | 978101285 |

Description

CLINEB PRO NEBULIZER - piston

High performing piston nebulizer for any kind of drug, ideal for intensive hospital and clinical use. Structure made of a high thermal and electrical insulated plastic in accordance with safety European standards. Long life dry piston compressor and bulb with a valve system allow to adjust the nebulization speed. Class IIa - Insulation class (EN60601-1): Class II Type B.

Made in Italy.

Technical Specifications

Operating Voltage: 230 V - 50 Hz
 Power consumption: 110 VA
 Max pressure: 300 kPa (3.0 bar)
 Air flow: 16 l/min
 Operating pressure: 95 kPa (0.95 bar)
 Operating air flow: 8.0 l/min
 Size - Weight: 23x19x25 cm - 2.2 kg
 Noise level: 60 dBA

Standard accessories

Bulb with valve
 Mouthpiece
 Nosepiece
 Adult and child masks with adapter
 Air tube
 Multilingual manual: GB, FR, IT, ES, DE