



Product Code	28065
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2835795
CND	Z12040210
EAN/UPC	6970384253290
PARAF	989370440

## Description

### MAGIC NEB WEARABLE MESH NEBULIZER

Pocket size and wearable ultrasonic nebulizer for treatment of respiratory diseases. Reliable quality, soundless operation and fine particle size. Self-cleaning mode: at the end of the treatment, the device mists pure water cleaning the mesh from medication residues.

- battery level indicator
- type-C charge port
- anti-corrosion material
- one-button operation

Supplied with adult mask, child mask, inhalation mouthpiece, built-in battery, type C cable and multilingual manual GB, FR, IT, ES, PT, DE, DK, RO, HR, GR.

## Technical Specifications

Ultrasonic frequency: approx. 110 kHz  
Nebulizing: rate = 0.2 ml/min  
Head life: 180 hours  
MMAD: < 5 µm  
FPF: = 60%  
Max./min. liquid volume: 6 ml / 2 ml  
Shutdown: automatic without liquid  
Power supply: built-in lithium battery 3.7 V, 300 mAh  
Operating time: > 60 min.  
Battery life: 500 cycles  
Size - weight: 54x36x44 mm - 45 g  
Liquid protection: IPX5