



Product Code	28066
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2952640
CND	Z12040210
EAN/UPC	8023279280661
PARAF	983277892

## Description

### MESH NEBULIZER

Small, handy, lightweight nebulizer for the treatment of respiratory diseases such as acute inflammation, tracheitis, bronchitis, swelling and pain in throat. It can operate with 2xAA standard batteries or USB cable.

- maximizes the benefits of medicine
- portable, handheld design (117 mm height)
- ultra-low noise, ultra-light only 104g
- one-button operation for easy use
- for adult and child: noiseless use and compact design allow to make inhalation even while children are sleeping

Supplied with adult and child mask, mouthpiece, 2xAA batteries, USB cable, carrying bag, multilingual manual and box (GB, IT, FR, ES, PT, DE, RO, GR).

## Technical Specifications

Power supply: 2xAA batteries or USB cable DC 5V  
 Ultrasonic frequency: approximately 110kHz  
 Nebulizer rate: =0.25ml / min  
 Max. / min. liquid volume: 10ml / 0.5ml  
 MMAD approximately 4µm  
 Nebulizing head life: 180 hours  
 Battery life: 60 minutes  
 Shutdown: automatic without liquid  
 Safety level: Type BF, Class II, IP 22