



Product Code 28070

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8054619101377

PARAF 970254660

## Description

### HUMI-RAINBOW HUMIDIFIER

A cold-mist ultrasonic humidifier with an attractive design and chromatic 7-led colours change

Mist output can be adjusted in direction (360° rotative head) and speed (adjustable knob)

The rear compartment can accommodate balsamic essences

## Technical Specifications

Power feeding: AC 110-240V / 50-60 Hz

Max mist output: 280 ml/h

Water Tank capacity: 2.8 l

Suitable for room size: min 10, max 25 m<sup>2</sup>

Max operating hours (non-stop): 10 hours

Size (individual box): Ø 19xh 30 cm - Weight: 1.4 kg