

COMP-A NEB PROFESSIONAL NEBULIZER 230V - 50 Hz



Product Code	28090
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011412
PARAF	977812243

Description

Piston compressor nebulizer for professional and home care use.
It embody all the dictates of modern aerosol therapy techniques.
Reliable and simple to use, also by inexperienced users.

Manual: GB, FR, IT, DE, ES.
Made in Italy.

Technical Specifications

- Operating voltage: 230 V - 50 Hz
- Size: 19.6 x 35.7 x h 18.5 cm
- Max pressure: 250 kPa
- Noise level: 55 dBA
- Operating pressure: 0-110 kPa
- Weight: 2.5 kg
- Air Flow: 15 l/min
- MDD class: IIa

Standard accessories

- Fasterjet nebulizer with valve system
- Nosepiece irrigation (polypropylene)
- Mouthpiece (polypropylene) with valve (silicone)
- Adult and child mask (PVC DEHP free)
- Air tube (PVC DEHP free) and spare filter
- Nasaljet