

"CORSIA PROFESSIONAL" NEBULIZER 230V - 50/60 Hz without manometer



Product Code	28092
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	237553
CND	Z12159002
EAN/UPC	8023279280920
GTIN	08023279280920
GMDN	35457

Description

CORSIA PROFESSIONAL NEBULIZER - piston - without manometer

Piston-type compressor nebulizers designed for continuous use. High performance with any type of drug, ideal for intensive hospital and clinic use. Recommended for atomizing antibiotics and bronchodilator drugs. Manufactured with high thermal and electric insulation plastic chassis in accordance with the latest European Safety regulations. The oil-free piston compressor has long durability and is equipped with the highly efficient HI-FLO jet nebulizer to guarantee quick and accurate drug delivery. Easy to handle and transport.

Made in Italy.

Technical Specifications

Operating Voltage: 220-230 V - 50/60 Hz - 184 VA
Max pressure: 250 kPa (2.5 bar)
Air flow: 16 l/min
Operating pressure: 130 kPa (1.30 bar)
Operating air flow: 5.2 l/min at 130 kPa
Nebulization rate: 0.50 ml/min with 4 ml solution
MMAD: 3.25 μ m
GSD: 3.45
Size: 255x190xh 105 mm
Weight: 2.2 kg
Noise level: approx 57 db

Standard accessories

HI-FLO Jet bulb
Mouthpiece
Nose piece
Adult and child mask
Air tube and air filter
Multilingual user manual and box: GB, FR, IT, ES, DE, PT, RO, RU*, PL*, GR*
*Box only