



Product Code	28094
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011474
PARAF	977812256
GMDN	12716

## Description

ATOMIZER NEBULIZER - 2 speeds

Piston-type compressor with 2 speeds for family use.  
With accessory compartment separated from cord compartment for improved hygiene.

Made in Italy.

## Technical Specifications

Operating voltage: 230 V - 50 Hz  
 Operating cycle: continuous  
 Operating pressure: speed 1: 45 kPa - speed 2: 60 kPa  
 Air Flow: speed 1: 8 l/min - speed 2: 11 l/min  
 Nebulization rate: speed 1: 0.2 ml/min - speed 2: 0.35 ml/min  
 Particle size: 80% < 5 µm  
 Max pressure: 1.85 bar / 185 kPa  
 Weight: 1.4 kg  
 Noise level: 56 db  
 Class of risk: IIa

## Standard accessories

Nebjet bulb with connector  
 Adult and child mask  
 Mouth piece and nose piece  
 Air tube and 3 air filter