



Product Code	28098
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2481965
CND	Z12159002
EAN/UPC	4015672112766

Description

OMRON C28P COMPRESSOR NEBULIZER - NE-C105-E

C28P nebulizer with Advanced Valve Technology has a breathenhanced valve system that ensures maximum drug flow during inhalation and no waste of medication during exhalation.

It can be used to treat chronic conditions or as relief from cough and cold, with medications or saline solution.

Designed for frequent use by family members and medical facilities.

- efficient treatment of lower airway conditions.
- long-life compressor.
- durable and easy to clean.

Technical Specifications

Capacity: min. 2 ml - max. 10 ml
 Operating mode: continuous use
 Particle size: approx. 2.65 µm
 Nebulizing rate: approx 0.12 ml/min
 Noise level: 55 dB
 Voltage: 230 V ~ 50 Hz, 140 VA
 Size - weight: 180x180x104 mm - 1.3 kg
 Safety: Type BF (applied part), Class II, IP21

Standard accessories

Nebulizer kit
 Child and adult mask
 PVC air tube (150 cm)
 Mouthpiece set
 5 spare air filters
 Power cable and adapter
 Carrying bag
 Multilingual manual GB, FR, IT, ES, DE, NL, TR, RU, Arabic.